


















Accessories

<p>Hanger for air mattress pump Attachment to hold air mattress pump. To be mounted on the footend of the bed.</p> <p>Part no.: PUMPBKT</p>	
<p>Care Assist Rail Designed to help patients safely in and out of bed.</p> <p>Part no.: 90702</p>	
<p>Push handle The push handle can be mounted on the foot end of the bed to improve the operation of the bed.</p> <p>Part no.: 90708 (Included in the product for FloorLine-i Plus, accessory for FloorLine-i)</p>	
<p>IV pole (adjustable in height) Adjustable IV pole that can be mounted on both ends of the bed.</p> <p>Part no.: IVP</p>	
<p>Extension kit Extension kit that extends the bed, with 175 mm to 2175 mm to accommodate more people.</p> <p>Part no.: 90709 (FloorLine LTC) Part no.: 90701 (FloorLine-i and FloorLine-i Plus)</p>	
<p>Oxygen Holder An attachment for placing an oxygen bottle in C size.</p> <p>Part no.: 90707</p>	
<p>Self help pole The self-help pole helps patients to get up in the bed and change positions.</p> <p>Part no.: SHP600 (FloorLine LTC) Part no.: 90700 (FloorLine i and FloorLine-i Plus)</p>	
<p>Side rails 1/2 side rail</p> <p>Part no.: 72791</p>	
<p>Side rails 1/4 side rail</p> <p>Part no.: 72790</p>	
<p>Side rails Split side rail (1/2 + 1/4 length side rail)</p> <p>Part no.: 72792</p>	
<p>Side rails 3/4 side rail</p> <p>Part no.: 72793</p>	
<p>Bumper The bumper protects the head of the bed and are mounted on the lower support bracket between the wheels.</p> <p>Part no.: 90705 (FloorLine LTC) Part no.: 90706 (FloorLine-i and FloorLine-i Plus)</p>	
<p>Swing-away trapeze The free-standing swing-away trapeze allows patient to reposition themselves. No tools needed for assembling.</p> <p>Part no.: 896</p>	
<p>Bed Raising Block (60 pcs/ box) The bed raising block is used for a simple way to raise or lower the seating or lying surface on beds, headboards, chairs or tables. May be stacked maximum 4 on top of each other.</p> <p>Part no.: 895</p>	
<p>Ring for Bed Raising Block The bed raising block can be complemented with the ring 895-1 for rounded and narrow legs (less than 35 mm) to obtain a solid filled block.</p> <p>Part no.: 895-1</p>	

” Bed solutions that provide reliability, high quality and innovation. Every step of the way.