



With the first ever cross-folding frame design, the neXus set the standard for what a rollator should be. The neXus is available in three colour choices and features a wide backstrap, ergonomic handles with optimal wrist positioning that optimizes the right amount of pronation and an ultimate comfort foam seat.

#### Cable-Free Braking System

The neXus is the only true-cable free rollator which not only means that there are no cables that can get caught on objects in and around the home to increase safety, but also ensures years of trouble-free use as cables cannot get loose within the frame itself. With the brakes in a locked position, the user can comfortably rest on our unique material foam seat.

#### Unique Frame Design

The neXus can stand independently, even in the folded position. The frame design allows the user ample space for an ergonomic walking pattern, or if needed, a slightly folded frame to pass through narrow doorways or

smaller areas. Because of this unique design, the stability is not compromised, ensuring the neXus is the most reliable and safe rollator on the market. The wide, double-band backstrap can be flipped to allow for two-position seating or used as a guide to assist with proper walking.

#### Lightweight and Easy to Fold

The cross-folding frame makes the neXus very easy to fold and store. The zippered soft bag comes as a standard option and can remain on the rollator, even in the folded position. With convenient accessories that can personalize the neXus to suit each individual need, the neXus is the only rollator you will ever need.

#### HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B  
117 43 Stockholm  
Phone: +46 8 665 35 00  
Fax: +46 8 665 35 10  
info.se@humancaregroup.com

#### HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street  
Port Melbourne, VIC 3207  
Phone: +61 1300 132 736  
Fax: + 61 (03) 9645 4954  
info.au@humancaregroup.com

#### HUMAN CARE CANADA

10-155 Colonnade Road  
Ottawa, ON K2E 7K1  
Phone: +613.723.6734  
Fax: +613.723.1058  
info.ca@humancaregroup.com

#### HUMAN CARE NEDERLAND

Elspeterweg 124  
8076 PA Vierhouten  
Phone: +31 577 412 171  
Fax: +31 577 412 170  
info.nl@humancaregroup.com

#### HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K  
Austin, TX 78754  
Phone: +512.476.7199  
Fax: +512.476.7190  
info.us@humancaregroup.com

## Accessories/Features



### Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 80702



### Seat handle

Included with rollator. Colour matches frame on blue and red (Black has red handle).



### Cup holder

Holds standard size bottle. Part no. 90394



### Slow-down brake

Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no. 4950



### Curb climber

Helps raise the rollator slightly to rest onto a curb. Included with rollator.



### Identification Tag

Attaches to rollator in seconds. Two design options. Part no. 90392 (red), 90393 (green)



### Cane holder

Simplifies transport of a cane or umbrella. Part no. 80703



### Flip backstrap

Extra wide, two-strap backstrap flips for two-position seating. Included with rollator.



### Zippered soft bag

Included with rollator. Maximum capacity 5kg/11lb. Part no. 93280



### Hard basket

Ideal for those great shopping adventures. Maximum capacity 4.5kg/10lb. Part no.: 4510



### Ergonomic handles

Cable-free handles feature an ergonomic grip and arthritic ball on handle for easy application.

This product is manufactured in accordance with all relevant ISO standards.

**FUNDED PRODUCT**

Ontario  
ADP (Assistive Devices  
Program) MW3D60045



Alberta Health Services  
AADL (Alberta Aids to Daily Living) C226

## Technical Specifications

| Description                        | Super Low                    | Low                          | Standard                               |
|------------------------------------|------------------------------|------------------------------|--|
| Part number                        | Black<br>Blue<br>Red         |                              |  |
|                                    | 470018<br>470018B<br>470018R | 470022<br>470022B<br>470022R | 470024<br>470024B<br>470024R           |
| Handle Height (Min - Max)          | 78cm - 85cm / 30.5" - 33.5"  | 78cm - 93cm / 30.5" - 36.5"  | 83cm - 100cm / 32.5" - 39.5"           |
| Seat Height <sup>1</sup>           | 48.3cm / 19"                 | 56cm / 22"                   | 61cm / 24"                             |
| Length                             | 66cm / 26"                   | 68.6cm / 27"                 | 70cm / 27.5"                           |
| Product Weight                     | 6.8kg / 15lb                 | 6.8kg / 15lb                 | 7.3kg / 16lb                           |
| <b>Specifications (All Models)</b> |                              |                              |  |
| Overall Width                      | 59.5cm / 23.4"               | Depth Folded                 | 25.4cm / 10"                           |
| Width between Handles              | 45cm / 17.5"                 | Wheel Size                   | 20cm x 3cm / 7.9" x 1.2"               |
| Seat Size                          | 40.6cm x 24.1cm / 16" x 9.5" | Distance between Rear Wheels | 47cm / 18.5"                           |
| Weight Capacity                    | 200kg / 441lb                | Bag Size (LxWxD when open)   | 40.6cm x 23cm x 17.8cm / 16" x 9" x 7" |

<sup>1</sup>Add 1"/2.5cm for seat pad

### Follow Us on Social Media



facebook.com/humancarecd



instagram.com/humancaregroup



linkedin.com/company/human-care-ab



pinterest.com/humancaregroup



youtube.com/channel/UC7fjINh3zTjRxOP1IDf2Da0