

Fellow Series

Fellow Supreme

The Fellow Supreme combines high quality ergonomic design with a low centre of gravity for increased safety and stability. Its vast handle height adjustment range and high weight capacity of 150kg make it easily suitable for most users.

With its unique fork construction, the Fellow Supreme is easy to maneuver. The soft braking system ensures that minimal hand pressure is required, allowing users with reduced hand functions to easily lock the wheels in place. The Fellow Supreme is also very easy to fold with a one-hand operation and can stand independently when folded.



Fellow Supreme

Fellow Classic

The Fellow Classic is a high quality rollator ideal for both indoor and outdoor use. Safe and reliable it is available in five models: standard, low, high, wide and narrow.

The Fellow Classic's unique fork construction allows for easy manoeuvrability. Like the Supreme, it features a soft braking system ideal for those users with limited hand mobility and an easy-to-use folding system that will accommodate a variety of users. Weight capacity 150kg.



Fellow Classic


Superior walking ergonomics

The forward facing handles allow the user to walk more upright which is easy on both the back as well as the shoulders. The Fellow series of rollators are also very easy to maneuver while also offering maximum stability. The sturdy seat flips up for added walking space or access to counters.
















Comfortable and easy to use

In order to offer you the best ergonomics, the handles are shaped to distribute the weight you put on them evenly. The soft braking system ensures that minimal hand pressure is required, allowing users with reduced hand functions to easily lock the wheels in place.

Technical Specifications								
Description	Part no.	Handle Height (min-max)	Seat height	Width	Width between handles	Wheel size	Weight	Max user weight
Fellow Supreme	85400	73 cm - 100 cm / 28.7" - 39.4"	59 cm / 23.2"	63 cm / 24.8"	44 cm / 17.3"	20 x 4 cm / 7.7"x1.7"	8.0 kg / 17.6 lbs	150 kg / 330.7 lbs
Fellow Classic std	88009	77 cm - 95 cm / 30.3" - 37.4"	59 cm / 23.2"	63 cm / 24.8"	44 cm / 17.3"	16 x 4 cm / 7.7"x1.7"	8.0 kg / 17.6 lbs	150 kg / 330.7 lbs
Fellow Classic low	88109	68 cm - 81 cm / 26.8" - 32.7"	51 cm / 19.9"	62 cm / 24.4"	44 cm / 17.3"		7.7 kg / 16.9 lbs	130 kg / 286.6 lbs
Fellow Classic wide	88309	77 cm - 95 cm / 30.3" - 37.4"	60 cm / 23.6"	69 cm / 27.2"	50 cm / 19.7"		8.7 kg / 19.2 lbs	150 kg / 330.7 lbs
Fellow Classic narrow	88509	77 cm - 95 cm / 30.3" - 37.4"	60 cm / 23.6"	60 cm / 23.4"	44 cm / 17.3"		8.1 kg / 17.8 lbs	130 kg / 286.6 lbs
Fellow Classic high	88609	82 cm - 100 cm / 32.3" - 39.4"	64 cm / 25.2"	63 cm / 24.6"	44 cm / 17.3"		8.1 kg / 17.8 lbs	130 kg / 286.6 lbs



Accessories

	Tray Practical for resting items. Part no, Fellow Classic: 08806 (Image 1)		Soft carry case ideal for personal items. Part no: 08808. (Cannot be used on Fellow Classic low)		Security Lock Flex cord lock that protects rollator from theft while unattended. Includes 2 keys. Part no: 91092
	Basket with handle Basket ideal for shopping. Part no, Fellow Supreme: 10667 (Image 2)		Back rest Gives support for the back while seated. Part no: 10668		Oxygen holder Part no, Fellow Classic: 08811. Part no, Fellow Supreme: 10671
	Basket with handle Basket ideal for shopping. Part no, Fellow Classic standard: 08805, Fellow Classic low: 88051, Fellow Classic narrow, wide and high: 88052. (Image 1)		Curb climber Helps raise the rollator onto a curb. Integrated in rollator for Fellow Supreme. Part no (Fellow Classic): 10661		Newspaper strap Can only be used on Fellow Supreme. Part no: 10679
	Slow down brake Prevents the rollator from moving too fast. Part no, right: 88121, Part no, left: 88122. Part no: 10625 (set)		IV-stand Can only be used with Fellow Classic. Part no: 08814		One-hand brake Allows control of brakes with one hand. Part no, Fellow Classic, right: 10677, left: 10678. Part no, Fellow Supreme, right: 10669, left: 10676
	Cane holder Simplifies the transportation of a cane or umbrella. Part no: 10603		Reflector Kit Includes 30 self-adhesive reflectors (10 white, 10 silver, 5 each non-flexible red and silver. Part no: 91090		Light Kit Battery operated, LED light mounts with either magnet or rubber fixture. Part no: 91091

Not available in all countries.